



## Speech-in-Noise Test For children, cochlear implant users and adults

### **BKB-SIN™** Speech-in-Noise Test

#### Developed for

- Children (ages 5 and up)
- Cochlear implant users
- Adults for whom the QuickSIN test is too difficult

#### Features

- Bamford-Kowal-Bench sentences in four-talker babble
- Uses easier sentences than the QuickSIN™
- Norms for children, adults and cochlear implant users

#### 18 equivalent List Pairs for estimating SNR Loss

- CD 1 (Standard BKB-SIN) has the target talker and background babble recorded on the same channel of the CD, at pre-recorded signal-to-noise ratios.
- CD 2 (Split Track BKB-SIN) has the target talker and background babble recorded on separate channels of the CD

#### Other uses

- Demonstrating the benefits of amplification
- Predicting performance with hearing aids in loud, noisy environments
- Assessing directional microphone performance
- Estimating children's performance for soft speech
- Screening for auditory processing disorders in children

#### Package includes

- BKB-SIN Test: Two-CD set
- Educational materials
- User manual
- Score sheets



PECIFICATION CAN BE SUBJECT TO CHANGE WITHOUT NOTICE

THE UK'S LEADING SOURCE OF AUDIOLOGICAL SOLUTIONS