

Kamplex K200D

Kamplex K200D Hearing Aid

The Kamplex K200D is a highly effective hearing aid designed to deliver high performance without undue complexity for a user and clinician.

Designed to be fitted as a standard BTE instrument or with a thin tube for high frequency losses and available in a choice of four colours.


Fitting using the K-Fit software requires minimal training to achieve a good result with simple toggles to fine-tune previous fittings.



Key Features

- 12 channel gain adjustment
- 12 channel layered noise reduction
- Up to 4 independent programs
- Data logging
- Earhook or thin tube fitting options
- 8 channel WDRC and MPO adjustment
- Reliant adaptive feedback cancellation
- Available in silver, charcoal grey, champagne beige or pearl white colour

Recommended Fitting Range



Accessories

Fitting Kit containing the below items:



- Ear Hooks
- Thin Tubes
- Open Fit Domes: S,M,L
- Occluded Single Dome: L
- Occluded Double Domes: S,M,L
- Sizing Tool

ZipTip (instant custom tip)

K-Fit Fitting Software

K-Fit is an all new fitting software, designed from the ground up, to ensure a fast fitting process. Featuring a logical fitting flow, K-Fit requires minimal training to quickly fit to the NAL NL-2 or DSL5 fitting formulae. Fine-tuning and adjustment is easy thereafter, with a series of toggles and on-screen displays.

Frequency Response



Reference Measurements

IEC 60118-7/ ANSI S3.22-2003 (2cc)	Earhook/ Thin Tube

IEC 60118-0 (Open Ear Simulator)	Earhook/ Thin Tube